

VETERINARY X-RAY

CMP-200



The CMP 200® High Frequency X-Ray Generator Series continues CPI's history of innovation by providing a high performance solution with enhanced features for the radiographic market that is unmatched in reliability.





動物用高頻X光機(32kw)



- ✓ CMP 200 高頻400kHz發生器系統
- ✓ 低幅射劑量，並取得高質量影像
- ✓ 內建犬、貓、禽類及特寵曝光條件
- ✓ 可自動變更大小焦點，或手動自行設定
- ✓ 兼容 14x17 或 17x17 英寸平板探測器
- ✓ kV、mA、Sec 可單獨調整

SOAR



Specification

Generator

Model/Maker	CMP 200 / CPI
Input Phase/Voltage	1Phase,220V/50A
Max output power	32kw
Generator type	Max.400kHz
kV _p Range	40-125kV _p
mA Range	10-400mA(10mA/step)
mAs Range	0.1-500mAs
Rotor Supply	Low Speed Starter
Exposure Timer Range	0.001s-6.3s
Anatomical programs (APR)	1024 techniques
Image Receptors	Up to 5
Technique Selection	1. kV/AEC option
	2. kV/mAs
	3. kV/mA/ms
	4. kV/cm thickness

Collimator

Max. Field Size	48 x 48cm at SID100cm
Auto Lamp Timer	Push Button, 30 sec.
Projection lamp	LED
Adjustment	Stepless

Tube

Model/Maker	E7239X / Canon
Focal spot	1.0/2.0
Max. Anode HU	100KJ(140KHU)
Anode target angle	16°

Tube Stand

Dimensions	2050(H)mm
Tube Rotation Angle	±90°
Vertical travel	±50cm

Table

Dimensions	1200x590x815(H)mm
Horizontal travel	±350mm



Careray
DR



Mikasa
X-Ray



JPI
X-Ray



Genoray
C-arm

總騰儀器股份有限公司

SOAR INSTRUMENT-TECH.CO.,LTD

新北市汐止區茄興街156號

TEL : 02-86486194 FAX : 02-86489678

E-MAIL : soar.jimmy@gmail.com

SOAR

CPI
Communications & Power Industries



Facebook



Official website