




medical imaging systems

WSR-40

High Frequency Inverter X-ray System



WSR-40




High Frequency Inverter X-ray System

Can reduce shooting time about 50% than general when WSR-40 system films human body by high frequency inverter method and because control console is small and light, is possible by wall tapestry and Radiographic is available by remote control switch conveniently that user films. And have A.P.R function basis.

SPECIFICATIONS

Component	Item	WSR-40
Console & Generator	kVp Range	40 ~ 125kV
	mA Range	50mA ~ 500mA
	mAs Range	0.5mAs ~ 630mAs
	Time Range	0.01 ~ 6.3 sec
	Self diagnosis	LED code display
	APR	576 kinds
	Output Frequency	60 KHz
X-Ray Tube (E7239X)	Tube Voltage	500mA/125kVp
	Focus	Small 1.0x1.0mm / Large 2.0x2.0mm
	Anode Heat	140,000HU
Collimator (CM-W1)	Ligting	More than 160Lux at SID 100cm
Tube Stand (TS-W1)	Ceiling Height	Approx. from 228cm above floor level
	U/D mov. Dis	Approx. 1370mm(380 to 1750mm)
	L/R mov. Dis	Approx. 1500mm Rail 2400mm
	Up / Down column rotation	± 180° with detent at each 90° position
	Horizontal rotation	± 180° with manual lock
Bucky Table 4 Way (BT-W1)	Table Dimensions	2000L x 728W x 710H
	Tabletop Sliding	Approx. more than ±375mm Longitudinal Approx. more than ±120mm Transverse
Bucky Stand (BS-W1)	U/D mov. Dis	Approx. 1500mm(320 to 1820mm)
H.T cable (HVC-6)	kVp Range	Max. 150kVp
Option Table	Fixed table (FBT-W3) 、 6 Way table (BT-W2) 、 Ceiling type (TC-W1)	



SOAR SOAR INSTRUMENT-TECH.CO.,LTD

No.156,Ln.20,Jiasing St.,Sijhih City,Taipei Jcounty 22172,Taiwan
Tel : 02-86486194 E-mail : soar.jimmy@gmail.com



Facebook



Official website